

De-identification process



USER USING
A BROWSER



AUTOMATIC DE-IDENTIFICATION BEFORE SENDING

Replaced items values to "ANONYMIZED"

- ✓ Patient's Name
- ✓ Patient's ID
- ✓ Patient's Birth Date
- ✓ Institution Name
- ✓ And more...

Deleted Items

- ✓ Physician ID Sequence
- ✓ Patient's Age
- ✓ Patient's Address
- ✓ And more...



https://



DE-IDENTIFICATION ON THE SERVER SIDE

PSEUDONYMIZATION FOLLOWING TABLE E.1-1 FROM PS 3.15 NEMA
DICOM STANDARD 2021B + CUSTOM RULES

Default Profile Kept Items (Original Value)

- ✓ Acquisition Date
- ✓ Patient's Sex
- ✓ Patient's Weight
- ✓ Study date
- ✓ Study Time

Common Custom Rules

- ✓ Study Description - Replaced by Study Trial Name
- ✓ Patient ID - Replaced by Patient's Study Code

Full Profile (No additional rule)



**GDPR
COMPLIANT**



IMAGING
MEDICAL EXPERT